


# BiOptic *Qexp*-MDx Direct SARS-CoV-2 Detection Kit-B

Direct multiplex one-step Reverse Transcription PCR testing for  
the detection of RNA from SARS-CoV-2



For Professional Use

 C903502-010, C903502-020, C903502-050,  
C903502-100 and C903502-200



BiOptic Inc.  
(23141) 4F., No.108-3, Minquan Rd., Xindian Dist.,  
New Taipei City, Taiwan (R.O.C.)



Medical Technology Promedt Consulting GmbH  
Altenhofstrasse 80, 66386 St. Ingbert, Germany

## CONTENTS

|   |           |
|---|-----------|
| <b>Product Information .....</b>                                | <b>3</b>  |
| Intended Use .....  | 3         |
| Product Description.....  | 4         |
| Contents and Storage .....                                      | 4         |
| Required materials not supplied.....                            | 5         |
| Assay limitations .....   | 6         |
| Warning and precautions.....                                    | 7         |
| Workflow .....  | 8         |
| <b>Preparation of Direct Lysis assay.....</b>                   | <b>9</b>  |
| Prepare Collection Buffer & RNA Lysis Buffer .....              | 9         |
| <b>Guidelines for buccal collection .....</b>                   | <b>10</b> |
| Buccal and Throat Swab (Brush).....                             | 10        |
| Buccal collection from the patient.....                         | 10        |
| <b>Preparation of one-step RT-PCR reaction .....</b>            | <b>11</b> |
| Prepare one-step RT-PCR Master Mix .....                        | 11        |
| Prepare buccal sample lysate.....                               | 11        |
| Prepare RT-PCR reaction .....                                   | 12        |
| <b>Perform RT-PCR using the thermal cycler .....</b>            | <b>12</b> |
| Setup and Run.....  | 12        |
| <b>Analysis and results .....</b>                               | <b>13</b> |
| The results of quality control (NTC and Positive Control) ..... | 13        |
| The result of Positive Control.....                             | 13        |
| Interpretation of the results .....                             | 14        |
| <b>Performance characteristics .....</b>                        | <b>15</b> |
| Limit of detection (LOD) .....                                  | 15        |
| Reactivity (Inclusivity) .....                                  | 16        |
| Cross-reactivity .....  | 17        |
| Interference .....  | 19        |
| Clinical evaluation .....                                       | 20        |
| <b>Documentation and support .....</b>                          | <b>21</b> |
| Related documentation.....                                      | 21        |
| Customer and technical support .....                            | 21        |
| <b>Symbol explanation .....</b>                                 | <b>22</b> |

## Product Information

### Intended Use

The BiOptic *Qexp*-MDx Direct SARS-CoV-2 Detection Kit-B includes the assays and positive control for the direct multiplex one-step reverse transcription PCR (RT-PCR) testing for the qualitative detection of RNA from SARS-CoV-2 in buccal specimens from individuals suspected of COVID-19 by their healthcare provider.

Results are for the identification of SARS-CoV-2 RNA. The SARS-CoV-2 RNA is generally detectable in the oral specimen during the acute phase of infection. Positive results are indicative of the presence of SARS-CoV-2 RNA; clinical correlation with patient history and other diagnostic information is necessary to determine patient infection status. Positive results do not rule out a bacterial infection or co-infection with other viruses. The agent detected may not be the definite cause of the disease. Laboratories within Europe and its territories are required to report all positive results to the appropriate public health authorities.

Negative results do not preclude SARS-CoV-2 infection and should not be used as the sole basis for patient management decisions. Negative results must be combined with clinical observations, patient history, and epidemiological information.

Testing with the BiOptic *Qexp*-MDx Direct SARS-CoV-2 Detection Kit-B is intended for use by qualified and trained clinical laboratory personnel specifically instructed and trained in the operation of the *Qamp<sub>mini</sub>*<sup>TM</sup> Thermal cycler, Bio-Fragment Analyzer (*Qsep*<sup>TM</sup> series) with its software (*Q-Analyzer*<sup>TM</sup>), and in vitro diagnostic procedures.

## Product Description

The Bioptic *Qexp*-MDx Direct SARS-CoV-2 Detection Kit-B includes the assays and positive control for the direct multiplex one-step RT-PCR testing for the qualitative detection of RNA from SARS-CoV-2 in buccal specimens from individuals suspected of COVID-19 by their healthcare provider. The Bioptic *Qexp*-MDx Direct SARS-CoV-2 Detection Kit-B includes the following components:

- **Direct Lysis assay**  
Sample lysis assay contains Collection Buffer, Proteinase K, RNase Inhibitor, and RNA Lysis Buffer.
- **Direct SARS-CoV-2 Kit assay**  
Direct SARS-CoV-2 multiplex one-step RT-PCR assay – multiplexed assays contain a primer set specific to the Envelope region (E) of the SARS-CoV-2 genome. Also, a primer set for internal control specific to the human RNase P gene.
- **Direct SARS-CoV-2 Kit Positive Control**  
The main component of positive control is synthetic RNA containing the target region of SARS-CoV-2 and the human RNase P.

## Contents and Storage

Table 1 Bioptic *Qexp*-MDx Direct SARS-CoV-2 Detection Kit-B

| Label   | Component           | Amount (μL) |        |        |         |         | Storage |
|---|---------------------|-------------|--------|--------|---------|---------|---------|
|   |                     | 10 rxn      | 20 rxn | 50 rxn | 100 rxn | 200 rxn |         |
| Direct Lysis assay                                |                     |             |        |        |         |         |         |
| CB  | Collection Buffer   | 1000        | 2000   | 5000   | 10mL    | 20mL    | 25°C    |
| PK  | Proteinase K        | 130         | 260    | 650    | 1300    | 2600    | -20°C   |
| RI  | RNase Inhibitor     | 2           | 4      | 10     | 20      | 40      | -20°C   |
| L   | RNA Lysis Buffer    | 200         | 400    | 1000   | 2000    | 4000    | -20°C   |
| Direct SARS-CoV-2 multiplex one-step RT-PCR assay |                     |             |        |        |         |         |         |
| 1   | Primer solution     | 50          | 100    | 250    | 500     | 1000    | -20°C   |
| 2   | Enzyme solution     | 25          | 50     | 125    | 250     | 500     | -20°C   |
| 3   | 2x RT-PCR Buffer    | 125         | 250    | 625    | 1250    | 2500    | -20°C   |
| PC  | Positive Control    | 8           | 16     | 40     | 80      | 160     | -20°C   |
| W   | Nuclease-free water | 100         | 200    | 500    | 1000    | 2000    | -20°C   |

- The Bioptic *Qexp*-MDx Direct SARS-CoV-2 Detection Kit-B is shipped on dry ice. The components of the kit should be arrived frozen except the Proteinase K, RNase Inhibitor and Enzyme solution (with glycerol).
- All components should be stored at -20°C immediately after unpackaging, except the Collection Buffer (store at 25°C).
- It is recommended to aliquot the reagents to avoid the risk of efficiency decrease by being thawed and frozen more than 2 times.

**Required materials not supplied**

Table 2 Direct sample collection tool (Not Provided)

| Item                             | Source        |
|----------------------------------|---------------|
| Direct sample collection tool    |               |
| Specimen Collection Brush (Swab) | any supplier* |

Table 3 Required Materials for BiOptic *Qexp*-MDx Direct SARS-CoV-2 Detection Kit-B (Not Provided)

| Item  | Source   |
|---|--|
| Consumables   |  |
| Disposable powder-free gloves and lab coats                                   | any supplier*  |
| Disposable pipette tips with the aerosol barrier                              | any supplier*  |
| 10% bleach or other appropriate cleaning solution that degrades nucleic acids | any supplier*  |
| RNaseOUT or 75% ethanol for bench cleaning                                    | any supplier*  |
| PCR tubes, plates, or strip tubes for the thermocycler being used             | any supplier*  |
| 1.5 mL microcentrifuge tubes, nuclease-free                                   | any supplier*  |
| Equipment   |  |
| Several micropipettes capable of pipetting volumes from 2µL to 1000µL         | any supplier*  |
| A cold block or ice   | any supplier*  |
| Vortex and centrifuge   | any supplier*  |
| Class II Biosafety cabinet, ideally in a BSL-2 containment facility           | any supplier*  |
| PCR workstation, for master mix preparing and setup                           | any supplier*  |
| <i>Qamp<sub>mini</sub></i> <sup>TM</sup> Thermal cycler (with program chips)  | BiOptic Cat. C310200   |
| Bio-Fragment Analyzers (instrument) ( <i>Qsep</i> series)                     | BiOptic Cat. C400100<br>BiOptic Cat. C100100<br>BiOptic Cat. C100001<br>BiOptic Cat. C100001-L |
| <i>Qsep</i> series cartridge kits   | BiOptic  |
| <i>Q-Analyzer</i> (software)  | BiOptic  |

\*The required materials not provided in this kit can be ordered from another major laboratory supplier.

### Assay limitations


- This kit is performed by using buccal mucosa samples. Nasopharyngeal samples could be detected as well, but have to be diluted 2 times and the LOD will be higher. Specimen types other than mucosa or nasopharyngeal swab could not be worked with this assay.
- The buccal mucosa samples with sputum may not be worked with this assay.
- The collection and storage of the samples must follow the appropriate procedures and conditions. Inappropriate uses may cause a decrease in the efficiency and sensitivity of the assay to detect the RNA sequence.
- False-negative results may cause by:
  - Improper sample collection or operation
  - Degradation of the virus RNA in the procedure
  - The presence of the RT-PCR inhibitors
- False-positive results may cause by:
  - Cross-contamination during patient samples collection or preparation
  - RNA contamination between patient samples
  - Cross-contamination by PCR amplicon
  - Contamination by improper preparation of Positive Control
- Negative and positive results are required to report to the appropriate public health authorities for proper patient management decisions.

### Warning and precautions

This Bioptic *Qexp*-MDx Direct SARS-CoV-2 Detection Kit-B is only used for *in vitro* diagnosis and should be operated on by medical professionals. The testing results should not be used for other purposes. To avoid cross-contamination and false-positive results, the Bioptic *Qexp*-MDx Direct SARS-CoV-2 Detection Kit-B workflow should use different areas or rooms. A dedicated PCR workstation for the preparation of the one-step RT-PCR assay would be ideal, and the analysis of the PCR product should be carried out in a separate area. The tube with Positive Control should be handled with great care to avoid contamination problems.

- Samples should always be treated as if infectious or biohazardous materials following safe laboratory procedures.
- Use personal protective equipment consistent with the guidelines for handling potentially infectious samples.
- Do not eat, drink, smoke, or apply cosmetic products in the work areas.
- Always use pipette tips with aerosol barriers. Tips should be sterile and DNases- and RNases- free.
- Reagents must be stored and handled as specified in Table 1.
- To prevent cross-contamination and false-positive results, be sure to use different pipettes for the preparation of one-step RT-PCR reagents and Master Mixes, for the preparation of Sample Lysates, and for pipetting the Positive Control.
- Modifications of assay reagents or assay protocols are not permitted.
- Do not use the kit after the expiration date.
- Dispose of waste in compliance with local, state, and federal regulations.
- After interpretation of results, all biological components should be discarded safely following the guideline revealed by the local authority.
- Always wash hands with soap and water or use hand sanitizer when complete.
- Laboratories and their territories are required to report all negative and positive results to the appropriate public health authorities.
- Positive results are indicative of the detection of SARS-CoV-2 RNA.

## Workflow



### ※ IMPORTANT

- To avoid false-positive results, always use different areas and pipettes for preparing reagents, adding positive control, adding buccal lysate, and analyzing PCR amplicon.
- Always use tips with the aerosol barrier to avoid contamination.
- Keep samples and reagents on the ice or freezing container during use.
- For preparing the reactions every time, include the following controls:
  - One positive control
  - One no template control (NTC)

The prepared RNA Lysis Buffer and Proteinase K can be stored at 4°C temporarily and should be used in 24hrs. It is recommended to prepare an appropriate amount of the above reagents and store them at 4°C, then do the test immediately or tomorrow. The reaction Master Mix should be prepared fresh on ice before use.

The RNA from buccal specimens could be released by using the required RNA Lysis Buffer. The RNA from buccal sample lysate is reverse transcribed into cDNA and amplified by the multiplex one-step RT-PCR on PCR thermal cycler. In the process, the primers anneal to a specific gene of SARS-CoV-2 and a human internal control as follows:

- E (envelope) protein
- Human RNase P

After RT-PCR, the PCR amplicon is analyzed by the *Qsep* series Bio-Fragment Analyzers, then provide the report automatically by the software *Q-Analyzer* (not required in the BiOptic *Qexp*-MDx Direct SARS-CoV-2 Detection Kit-B).



## Preparation of Direct Lysis assay

### Prepare Collection Buffer & RNA Lysis Buffer

1. Thaw RNA Lysis Buffer and Nuclease-free Water on ice.
2. Take out the Proteinase K and RNase Inhibitor and keep them cold in a freezing container.
3. Centrifuge briefly to collect the liquid at the bottom of the tube.
4. Prepare the required amount of clean 1.5 mL microcentrifuge tube.
5. Aliquot 13  $\mu\text{L}$  Proteinase K in a 1.5 mL microcentrifuge tubes for every sample. (e.g., prepare 5 tubes for 5 buccal samples, and aliquot 13  $\mu\text{L}$  Proteinase K in every tube.)
6. Add 100  $\mu\text{L}$  Collection Buffer in every 1.5 mL microcentrifuge tube containing 13  $\mu\text{L}$  Proteinase K.
7. Mix well by pulse vortex, then centrifuge briefly to collect the liquid at the bottom of the tube.
- ※ If there is precipitation in Collection Buffer, please keep it at room temperature until it is completely dissolved before use.

Table 4 Preparation of the RNA Lysis Buffer

| Component        | Volume per sample  | Volume for n samples        | Volume for 94 samples |
|------------------|--------------------|-----------------------------|-----------------------|
| RNase Inhibitor  | 0.2 $\mu\text{L}$  | 0.2 $\mu\text{L} \times n$  | 18.8 $\mu\text{L}$    |
| RNA Lysis Buffer | 20 $\mu\text{L}$   | 20 $\mu\text{L} \times n$   | 1880 $\mu\text{L}$    |
| TOTAL            | 20.2 $\mu\text{L}$ | 20.2 $\mu\text{L} \times n$ | 1898.8 $\mu\text{L}$  |

8. Gently mix the RNA Lysis Buffer by inverting or vortex, then centrifuge briefly to collect the liquid at the bottom of the tube.
9. The aliquoted Proteinase K and prepared RNA Lysis Buffer can be stored at 4°C and should be run out in 24hrs.

## Guidelines for buccal collection

### IMPORTANT


- Before buccal collection, rinse the mouth with clean water for 10 seconds to remove microorganisms, and food and drink remnants.
- Ensure that there was no eating, drinking, smoking, chewing gum, brushing teeth, or use of mouthwash for at least 1 hour after rinsing the mouth.
- Before specimen collection, cough slightly for 5 times to increase the detectable viral load but avoid coughing up sputum.

### Buccal and Throat Swab (Brush)


Flocked, foam, or nylon swab with plastic shafts and the transport tube. Do not use calcium alginate swabs or swabs with wooden shafts, they may contain substances which could impair the test.

### Buccal collection from the patient

1. Fill in the form on the collection tube.
2. Open the Specimen Collection Brush (Swab).
3. Rub and rotate the brush/swab on the buccal mucosa and gum 10 times on each side to collect the virus and cells.
4. Put the brush/swab back to the collection tube. (Transportation of the COVID-19 samples please follow the Guidelines released by CDC)



5. Snap the brush/swab (optional), then insert it into the Collection Buffer & Proteinase K mixture to release the virus and cells.



6. Ensure that the buccal sample is mixed well by stirring or vortex, then throw out the brush/swab and centrifuge briefly to collect the liquid at the bottom of the tube.
7. Incubate at room temperature for 1 min.  
(The buccal mixture could be stored at 4°C for 24 hours.)

### Preparation of one-step RT-PCR reaction

#### Prepare one-step RT-PCR Master Mix

1. Thaw the Primer Solution and 2x RT-PCR buffer on ice.
2. Gently vortex the Primer Solution and 2x RT-PCR buffer, then centrifuge briefly to collect the liquid at the bottom of the tube.
3. Prepare the reaction Master Mix in the PCR tube or 1.5 mL microcentrifuge tube:

Table 5 Preparation of RT-PCR reagents

| Component        | Volume per sample | Volume for n samples plus 2 controls | Volume for 94 samples plus 2 controls |
|------------------|-------------------|--------------------------------------|---------------------------------------|
| Primer solution  | 5 µL              | 5 µL x (n+2)                         | 480 µL                                |
| Enzyme solution  | 2.5 µL            | 2.5 µL x (n+2)                       | 240 µL                                |
| 2x RT-PCR buffer | 12.5 µL           | 12.5 µL x (n+2)                      | 1200 µL                               |
| TOTAL            | 20 µL             | 20 µL x (n + 2)                      | 1920 µL                               |

4. Mix by pipetting or gently invert the reaction Master Mix, then centrifuge briefly to collect the liquid at the bottom of the tube.
5. Transfer every 20 µL Master Mix in a PCR tube.
- ※ The reaction Master Mix should be prepared fresh on ice before use.

#### Prepare buccal sample lysate

1. Transfer 20 µL buccal mixture into a clean PCR tube.
2. Heated to 95°C for 5 mins to inactivate Proteinase K, then cooled down to 4°C by *Qamp<sub>mini</sub>*<sup>TM</sup> with program chip A.
3. Add 20 µL RNA Lysis Buffer in a buccal mixture immediately and mix well by pulse vortex, then centrifuge briefly to collect the liquid at the bottom of the tube and incubate for 10 mins at room temperature.

### Prepare RT-PCR reaction

Table 6 The preparation list of all RT-PCR reaction

| Component                  | Volume per reaction |                 |                  |
|----------------------------|---------------------|-----------------|------------------|
|                            | NTC                 | Sample reaction | Positive control |
| one-step RT-PCR Master Mix | 20 $\mu$ L          | 20 $\mu$ L      | 20 $\mu$ L       |
| Buccal sample lysate       | -                   | 5 $\mu$ L       | -                |
| Nuclease-free water        | 5 $\mu$ L           | -               | -                |
| Positive control           | -                   | -               | 5 $\mu$ L        |
| TOTAL                      | 25 $\mu$ L          | 25 $\mu$ L      | 25 $\mu$ L       |

- ※ Do not use the same pipette to add water, buccal sample lysates, and Positive Controls. Do not add Positive Control in the same area as reagents preparing. Avoid producing lots of bubbles in the PCR tube.
  - ※ Always prepare the NTC first, and the Positive Control the last to avoid cross-contamination.
1. For NTC, transfer 5  $\mu$ L nuclease-free water in the prepared Master Mix.
  2. For sample reaction, transfer 5  $\mu$ L buccal sample lysate in the prepared Master Mix.
  3. For Positive Control, transfer 5  $\mu$ L Positive Control in the prepared Master Mix after thawing it on ice.
  4. Mix by invert briefly or vortex gently after addition, then centrifuge briefly to collect the liquid at the bottom and remove the bubble if there it was.

### Perform RT-PCR using *Qamp<sub>mini</sub>*<sup>TM</sup>

Before using the thermal cycler, portable *Qamp<sub>mini</sub>*<sup>TM</sup> Thermal cycler, please read the provided operation manual of *Qamp<sub>mini</sub>*<sup>TM</sup> Thermal cycler and pay attention to the safety information. To guarantee problem-free operation, please follow the instructions and safety precautions to ensure safe operation.

### Setup and Run

1. The thermal protocol in program chip B as below:

Table 7 The program of the Direct SARS-CoV-2 multiplex one-step RT-PCR assay

| Step                  | Number of cycles | Temperature ( $^{\circ}$ C) | Time       |
|-----------------------|------------------|-----------------------------|------------|
| Reverse transcription | 1                | 50                          | 15 minutes |
| Activation            | 1                | 95                          | 10 minutes |
| Denaturation          | 44               | 95                          | 15 seconds |
| Anneal                |                  | 60                          | 20 seconds |
| Extension             |                  | 72                          | 20 seconds |
| Extension             | 1                | 72                          | 5 minutes  |
|                       | 1                | 4                           | hold       |

2. After the RT-PCR is complete, transfer the PCR tube to the CGE system for analysis.

## Analysis and results

For better sensitivity and resolution, please use capillary gel electrophoresis (CGE) system, such as BiOptic *Qsep* series Bio-Fragment Analyzers, to analyze the RT-PCR amplicon. Before using the capillary gel electrophoresis system, please read the operation manual which the manufacturer provided and pay attention to the safety information. For detailed information on the cartridge and BiOptic's Capillary Gel Electrophoresis system, please see the provided operation manual or visit our website for more information.

Table 8 The size of the target sequences


|             | E      | RNase P |
|-------------|--------|---------|
| Target size | 129 bp | 170 bp  |

### The results of quality control (NTC and Positive Control)

- Ensure that there is no signal of the target sequences revealed in the NTC.
  - If the NTC has any signal, the reagents must be contaminated in the procedure.
- Ensure that all signals of the target sequences are shown in the Positive Control.
  - If the Positive Control did not show any signal, it may cause by the degradation of the reagents. Please contact the distributor in your area.

### The result of Positive Control

The result shown below was run by BiOptic's Capillary Gel Electrophoresis system, Bio-fragment analyzer (*Qsep*<sub>100</sub><sup>TM</sup>) with BiOptic's standard cartridge (S2).



### Interpretation of the results

The interpretation could be performed automatically by the software *Q-Analyzer* for the *Qsep* series (Table 9).

Table 9 Result interpretation for buccal samples

| E gene | RNase P | Result                             | Action  |
|--------|---------|------------------------------------|---|
| —      | —       | N/A                                | Repeat test. If the repeat remains invalid, consider a new specimen.  |
| —      | +       | SARS-CoV-2<br>NOT Detected         | Report results to the healthcare provider. Consider testing for other viruses.  |
| +      | —       | SARS-CoV-2<br>Presumptive Positive | Repeat test. If the repeat result remains inconclusive, additional confirmation testing should be conducted if clinically indicated. Report results to the healthcare provider and appropriate public health authorities. |
| +      | +       | Positive SARS-CoV-2                | Report results to the healthcare provider and appropriate public health authorities.  |

## Performance characteristics

The performance of the BiOptic *Qexp*-MDx Direct SARS-CoV-2 Detection Kit-B was evaluated by determining the limit of detection (LOD), cross-reactivity, and characterizing the impact of interfering substances as follows.

### Limit of detection (LOD)

The LOD study revealed the lowest SARS-CoV-2 viral concentration that can be detected by the BiOptic *Qexp*-MDx Direct SARS-CoV-2 Detection Kit-B in the buccal specimen. The negative buccal specimen was spiked with AccuPlex SARS-CoV-2 non-replicative recombinant viruses, then follow by the BiOptic *Qexp*-MDx Direct SARS-CoV-2 Detection Kit-B workflow.

Preliminary LOD was performed first and tested at four different concentrations (100, 50, 20, 10 copies in a reaction) in quintuplicate (Table 10). Then, the preliminary LOD was confirmed by testing 20 replicates at the concentration of 50 and 20 copies/reaction including quality control (negative samples and NTC in triplicate). The LOD ultimately was determined to be 20 copies in a reaction (Table 11).

Table 10 The preliminary LOD of four different concentrations

| Effective concentration<br>(copies/reaction) | Detection Rate |            |
|--|----------------|------------|
|  | E              | RNase P    |
| 100  | (5/5) 100%     | (5/5) 100% |
| 50   | (5/5) 100%     | (5/5) 100% |
| 20   | (5/5) 100%     | (5/5) 100% |
| 10   | (2/5) 40%      | (5/5) 100% |

Table 11 The 20 replicates LOD testing at 50, 20 copies in a reaction

| Effective concentration<br>(copies/reaction) | Detection Rate |              |
|--|----------------|--------------|
|  | E              | RNase P      |
| 50   | (20/20) 100%   | (20/20) 100% |
| 20   | (20/20) 100%   | (20/20) 100% |
| Negative                                     | (0/3) 0%       | (3/3) 100%   |
| NTC  | (0/3) 0%       | (0/3) 0%     |

### Reactivity (Inclusivity)

The designed primer sequences for detecting the SARS-CoV-2 E gene were aligned to 11 SARS-CoV-2 genomes from the database of the NCBI. The primers for the SARS-CoV-2 E gene have 100% identity to all 11 strains.

Table 12 The identity for primer pairs with SARS-CoV-2 sequences

| Sequence ID | Strain   | Identity for E target (%) |                 |
|-------------|--|---------------------------|-----------------|
| MZ576185.1  | Isolate SARS-CoV-2/human/USA/LA                            | 100                       | can be detected |
| MZ574958.1  | Isolate SARS-CoV-2/human/USA/TX                            | 100                       | can be detected |
| MZ573385.1  | Isolate SARS-CoV-2/human/UZB/UZB                           | 100                       | can be detected |
| NC_045512.2 | Isolate SARS-CoV-2/Wuhan-Hu-1                              | 100                       | can be detected |
| OU407641.1  | Isolate hCoV-19/Switzerland/ZH                             | 100                       | can be detected |
| MZ544374.1  | isolate SARS-CoV-2/human/KOR/KD-SARS2-2/2020               | 100                       | can be detected |
| MZ277385.1  | Isolate SARS-CoV-2/human/TWN/CGMH                          | 100                       | can be detected |
| OU424957.1  | SARS-CoV-2 genome assembly/UK, chromosome: 1               | 100                       | can be detected |
| OU398468.1  | SARS-CoV-2 genome assembly/Germany, chromosome: 1          | 100                       | can be detected |
| LC593804.1  | SARS-CoV-2 TKYE622388_2020 genomic RNA/JP, complete genome | 100                       | can be detected |
| MZ572206.1  | Isolate SARS-CoV-2/human/IND/CWS_0091_54742087/2021        | 100                       | can be detected |



### Cross-reactivity

*In silico* analysis of the following organisms

Table 13 Organisms used for *in silico* cross-reactivity analysis

| High priority pathogen                               |   |
|--|---|
| Human coronavirus 229E                               | Human coronavirus NL63  |
| Human coronavirus OC43                               | SARS-coronavirus  |
| Human coronavirus HKU1                               | MERS-coronavirus  |
| High priority organisms likely                       |   |
| Human Metapneumovirus (hMPV)                         | <i>Streptococcus pneumoniae</i>   |
| Parainfluenza virus 1-4                              | <i>Streptococcus pyogenes</i>   |
| Influenza A & B                                      | <i>Bordetella pertussis</i>   |
| <i>Enterovirus</i> (e.g. EV68)                       | <i>Mycoplasma pneumoniae</i>  |
| Respiratory syncytial virus                          | <i>Pneumocystis jirovecii</i> (PJP)   |
| Rhinovirus   | <i>Candida albicans</i>   |
| Adenovirus type 1, 7                                 | <i>Cytomegalovirus</i>  |
| Epstein Barr Virus                                   | <i>Measles morbillivirus</i>  |
| Mumps Virus genotype G                               | <i>Corynebacterium sp.</i> WS2071   |
| <i>Escherichia coli</i>                              | <i>Lactobacillus sp.</i>  |
| <i>Moraxella catarrhalis</i>                         | <i>Neisseria meningitidis</i>   |
| <i>Neisseria sp.</i>                                 | Human <i>Metapneumovirus</i> (hMPV)   |
| <i>Chlamydia pneumoniae</i>                          | <i>Pseudomonas aeruginosa</i>   |
| <i>Haemophilus influenzae</i>                        | <i>Staphylococcus epidermis</i>   |
| <i>Legionella pneumophila</i>                        | <i>Streptococcus salivarius</i>   |
| <i>Mycobacterium tuberculosis</i>                    | <i>Streptococcus pneumoniae</i>   |
| <i>Staphylococcus aureus</i><br>(Protein A producer) | Pooled human nasal wash - to represent diverse microbial flora in the human respiratory tract |

Among the organisms in Table 13, blast analysis showed  $\leq 80\%$  similarity to the E primers, except SARS-coronavirus. The analysis of SARS-coronavirus showed high similarity to the forward and reverse primers for the E gene. So, we choose the position where there is a mutation at the 3' end of the E gene forward primer of SARS-CoV-2. In addition, there is no case about SARS coronavirus now worldwide. Therefore, it will not affect the analysis and interpretation of the results.

*In vitro* analysis of the following organisms

Negative specimen sample spiked in NATrol™ Respiratory Verification Panel as a mimic specimen for Direct SARS-CoV-2 RT-PCR test. PC means provided SARS-CoV-2 Positive Control and NTC means no template control.

Table 14 The results of the cross-reactivity test

| Pathogen Name              | Target Signal |    | Quality Control |     | SARS-CoV-2   |
|----------------------------|---------------|----|-----------------|-----|--------------|
|                            | E             | IC | PC              | NTC |              |
| Influenza AH1              | -             | +  | +               | -   | not detected |
| Influenza AH3              | -             | +  | +               | -   | not detected |
| Influenza A H1N1pdm**      | -             | +  | +               | -   | not detected |
| Influenza B                | -             | +  | +               | -   | not detected |
| <i>Metapneumovirus</i> *** | -             | +  | +               | -   | not detected |
| RSV A                      | -             | +  | +               | -   | not detected |
| Rhinovirus                 | -             | +  | +               | -   | not detected |
| Parainfluenza 1            | -             | +  | +               | -   | not detected |
| Parainfluenza 2            | -             | +  | +               | -   | not detected |
| Parainfluenza 3            | -             | +  | +               | -   | not detected |
| Parainfluenza 4            | -             | +  | +               | -   | not detected |
| Adenovirus 3               | -             | +  | +               | -   | not detected |
| Coronavirus NL63           | -             | +  | +               | -   | not detected |
| Coronavirus 229E           | -             | +  | +               | -   | not detected |
| Coronavirus OC43           | -             | +  | +               | -   | not detected |
| Coronavirus HKU-1          | -             | +  | +               | -   | not detected |
| <i>M. pneumoniae</i>       | -             | +  | +               | -   | not detected |
| <i>C. pneumoniae</i>       | -             | +  | +               | -   | not detected |
| <i>B. pertussis</i>        | -             | +  | +               | -   | not detected |

\*\* Please note that although similar in nomenclature, this is a 2009 H1N1 pandemic Influenza strain and does NOT correlate with the seasonal 2009 Influenza strains found in the Fludb.org database. For reference, the NCBI Taxon IDs for the seasonal Influenza strains listed in the Fludb.org database are A/New York/01/2009 (H1N1) -666252; B/New York/01/2009 -664512; A/New York/02/2009 (H1N1)- 666298; and A/New York/03/2009 (H3N2) -659637.

\*\*\* "This product is sold by Zeptomatrix Corporation under license from Vironovative B. V under patent applications, including U.S. Patent Applications 10/371,099 and 10/371,12 and any patents that issue from applications related to PCT/NL02/00040 and PCT/US03/05271."

### Interference

The test is to verify the influence of potentially interfering substances on the BiOptic *Qexp*-MDx Direct SARS-CoV-2 Detection Kit-B. The material is prepared by a mixture with the potentially interfering substance in the buccal samples. Buccal samples without interfering substances were included as experimental sample controls, and all samples were prepared by spiking AccuPlex™ SARS-CoV-2 Verification Panel in the 2.5 x LOD. The analysis was shown in the following Table 15.

To minimize the effect of the potential inhibitor from specimens, the patient should rinse the mouth with clean water for 10 seconds first to remove microorganisms, and food and drink remnants. And the patient should not eat, drink, smoke, chew gum, brush teeth, or use mouthwash for at least 1 hour after rinsing the mouth before specimen collection. The collection buffer and RNA lysis buffer also have the effect minimize the interference of the inhibitor.

Table 15 Analysis results of potential interference substances

| Interference Substances                                 | Components                   | Concentration | Result       |
|---|------------------------------|---------------|--------------|
| NAZAL SPRAY   | Benzalkonium Chloride        | 1 µg/mL       | 3/3 Positive |
|   | Chlorpheniramine Maleate     | 50 µg/mL      | 3/3 Positive |
|   | Naphazoline Hydrochloride    | 5 µg/mL       | 3/3 Positive |
| Otrivin Menthol 0.1%, nasal drops in Metered-Dose-Spray | Xylometazoline hydrochloride | 10 µg/mL      | 3/3 Positive |
| Otrivin Anti-Allergy Nasal Spray                        | Fluticasone propionate       | 5 µg/mL       | 3/3 Positive |
| Bisoucan Nasal Spray "C.M."                             | Diphenhydramine HCl          | 20 µg/mL      | 3/3 Positive |
|   | Naphazoline HCl              | 5 µg/mL       | 3/3 Positive |
|   | Procaine HCl                 | 100 µg/mL     | 3/3 Positive |
| Mucin from porcine stomach                              |                              | 2.5 mg/mL     | 3/3 Positive |
| NICORETTE® Nicotine Inhaler                             | Nicotine                     | 0.03 mg/mL    | 3/3 Positive |
| Fucole Anti-Inflammatory Spray                          | Benzydamine Hydrochloride    | 7.5 µg/mL     | 3/3 Positive |
| Colgate® Toothpaste                                     |                              | 0.5% v/v      | 3/3 Positive |
| LISTERINE® Healthy White Mouthwash                      |                              | 5% v/v        | 3/3 Positive |
| Whole blood   |                              | 1% v/v        | 3/3 Positive |

### Clinical evaluation

The clinical evaluation was performed by comparing the EUA authorized molecular real-time RT-PCR assay (TBG ExProbe™ SARS-CoV-2 Testing Kit, Cat. No. 68020) with BiOptic *Qexp*-MDx Direct SARS-CoV-2 Detection Kit-B.

The samples for real-time RT-PCR assay were contrived nasopharyngeal swab specimens, spiked with AccuPlex SARS-CoV-2 non-replicative recombinant viruses as the positive specimen (1.5 x LOD); Invitrogen DEPC-treated water as the negative specimen, then all were extracted by QIAamp Viral RNA Mini Kit. The real-time RT-PCR assay was performed, and the results were analyzed by Roche LightCycler 480 II as below.

The samples for BiOptic Direct RT-PCR assay were contrived buccal swab specimens, spiked with the same AccuPlex recombinant viruses as the positive specimen (1.5 x LOD); Invitrogen DEPC-treated water as the negative specimen. Performing Direct RT-PCR by TurboCycler and the results were analyzed by BiOptic *Qsep* series as below.

Table 16 Clinical evaluation using contrived buccal & nasal specimens

| Buccal samples \ Nasal samples                                   |                       | ExProbe™ SARS-CoV-2 Testing Kit |          |       |
|--|-----------------------|---------------------------------|----------|-------|
|  |                       | Positive <sup>a</sup>           | Negative | Total |
| BiOptic <i>Qexp</i> -MDx<br>Direct SARS-CoV-2<br>Detection Kit-B | Positive <sup>a</sup> | 30                              | 0        | 30    |
|  | Negative              | 0                               | 30       | 30    |
|  | Total                 | 30                              | 30       | 60    |
| Positive Percent Agreement                                       |                       | 100% (88.43-100%) <sup>b</sup>  |          |       |
| Negative Percent Agreement                                       |                       | 100% (88.43-100%) <sup>c</sup>  |          |       |

<sup>a</sup> 1.5 x LOD is about 30 copies in a reaction.

<sup>b</sup> 95% CI of the Sensitivity

<sup>c</sup> 95% CI of the Specificity

Table 17 Clinical evaluation using contrived nasal specimens

| Nasopharyngeal swab  |                       | ExProbe™ SARS-CoV-2 Testing Kit |          |       |
|--|-----------------------|---------------------------------|----------|-------|
|  |                       | Positive <sup>a</sup>           | Negative | Total |
| BiOptic <i>Qexp</i> -MDx<br>Direct SARS-CoV-2<br>Detection Kit-B | Positive <sup>a</sup> | 30                              | 0        | 30    |
|  | Negative              | 0                               | 30       | 30    |
|  | Total                 | 30                              | 30       | 60    |
| Positive Percent Agreement                                       |                       | 100% (88.43-100%) <sup>b</sup>  |          |       |
| Negative Percent Agreement                                       |                       | 100% (88.43-100%) <sup>c</sup>  |          |       |

<sup>a</sup> 1.5 x LOD is about 30 copies in a reaction.

<sup>b</sup> 95% CI of the Sensitivity

<sup>c</sup> 95% CI of the Specificity

## Documentation and support

### Related documentation

Table 18 The related documents of BiOptic's kits or instruments












| Document   | Link  |
|--|---|
| <i>Qamp<sub>mini</sub></i> <sup>TM</sup> Thermal cycler Operation Manual                         | <a href="https://apps.bioptic.com.tw/webdl/Instrument/F0026-Qamp%20mini%20PCR%20Manual%20ENG-B.pdf">https://apps.bioptic.com.tw/webdl/Instrument/F0026-Qamp%20mini%20PCR%20Manual%20ENG-B.pdf</a>   |
| <i>Qsep</i> Series Software Installation Quick Start   | <a href="https://apps.bioptic.com.tw/webdl/Instrument/F0002-Software%20Installation%20Quick%20Start%20-ENG-E.pdf">https://apps.bioptic.com.tw/webdl/Instrument/F0002-Software%20Installation%20Quick%20Start%20-ENG-E.pdf</a>                   |
| <i>Qsep<sub>1</sub></i> <sup>TM</sup> & <i>Qsep<sub>100</sub></i> <sup>TM</sup> Operation Manual | <a href="https://apps.bioptic.com.tw/webdl/Instrument/F0044_Qsep1%20and%20Qsep100%20Operation%20Manual-Hardware-ENG-C.pdf">https://apps.bioptic.com.tw/webdl/Instrument/F0044_Qsep1%20and%20Qsep100%20Operation%20Manual-Hardware-ENG-C.pdf</a> |
| <i>Qsep<sub>1</sub>-Lite</i> <sup>TM</sup> Operation Manual                                      | <a href="https://apps.bioptic.com.tw/webdl/Instrument/F0052_Qsep1-Lite%20Operation%20Manual-Hardware-ENG-B.pdf">https://apps.bioptic.com.tw/webdl/Instrument/F0052_Qsep1-Lite%20Operation%20Manual-Hardware-ENG-B.pdf</a>                       |
| <i>Qsep<sub>400</sub></i> <sup>TM</sup> Operation Manual   | <a href="https://apps.bioptic.com.tw/webdl/Instrument/F0043_Qsep400%20Operation%20Manual-%20Hardware%20-ENG-E.pdf">https://apps.bioptic.com.tw/webdl/Instrument/F0043_Qsep400%20Operation%20Manual-%20Hardware%20-ENG-E.pdf</a>                 |
| 2020 BiOptic Brochure  | <a href="http://apps.bioptic.com.tw/webdl/Catalogue%20&amp;%20Flyers/2020_BiOptic%20catalog_low_resolution.pdf">http://apps.bioptic.com.tw/webdl/Catalogue%20&amp;%20Flyers/2020_BiOptic%20catalog_low_resolution.pdf</a>                       |

### Customer and technical support

For more information about the detection kit or instruments, please visit BiOptic's website: <https://www.bioptic.com.tw/?lang=EN>.

## Symbol explanation

Table 19 Explanation of the symbols used in the labeling

| Symbol  | Symbol Title  | Explanatory Text   |
|---|---|--|
|    | CE marking Conformité Européene Notified Body Reference no. ### | The requirements for accreditation and market surveillance relation to the marketing of products |
|    | Manufacturer  | Indicates the medical device manufacturer  |
|    | Date of manufacture   | Indicates the date when the medical device was manufactured                                      |
|    | Authorized Representative in the European Community             | Indicates the authorized representative in the European Community                                |
|   | In vitro diagnostic medical device                              | Indicates a medical device that is intended to be used as an in vitro diagnostic medical device  |
|  | Catalog number  | Indicate the manufacturer' s catalog number to identify the medical device                       |
|  | Batch code  | Indicates the manufacturer' s batch code to identify the batch or lot                            |
|  | Use by  | Indicates the date after which the medical device is not to be used                              |
|  | Upper limit of temperature                                      | Indicates the upper limit of temperature to which the medical device can be safely exposed       |
|  | Consult instructions for use                                    | Indicates the need for the user to consult the instructions for use                              |
|  | Keep away from sunlight   | Indicates a medical device that needs protection from light sources                              |