

AccuGO™ Pfu


Proofreading DNA polymerase




AccuGO pfu PCR Kit

is a high quality, high performance PCR proofreading kit. The new generation of Pfu enzyme is engineered to provide high efficiency, high accuracy and high yield PCR products for widely use in cloning, gene expression and sequencing.....etc.. The AccuGO PCR Kit provides both short extension time and high specificity. High GC content template or fragment up to 10kb can also be produced with high yield.

(a) High specificity and high yield



(b) Successful PCR with >70% GC containing template



(c) High yield amplification of long PCR products



- Specialized engineered Pfu enzyme for **high fidelity products**
- Short extension time : **15-30s/ kb**
- **High yield PCR products**: Compare with standard Taq kit
- High yield for long fragment PCR product: **up to 10kb**
- Good performance in **GC rich** template
- PCR product with **blunt ends**